

VATS Trainer for Explants

LEFT-SIDED



**Bringing
Practice to the
Practice of
Medicine®**

Provides a realistic human thoracic context for lifelike VATS training with porcine tissue explants

Derived from patient CT data, the VATS Trainer for Explants Left-Sided is a hemi-thorax with all pertinent landmarks for thoracoscopic lung surgeries. The trainer presents a patient in a lateral decubitus position with exposed anterior and posterior thoracic wall. The replaceable skin/muscle element supports trocars and is comprised of our renowned mimetic tissue. Ribs may be palpated beneath the soft tissue for procedural port placement. The trainer is equipped with a tray to capture tissue effluence and is easily disassembled for cleaning.

TO ORDER CALL: +1.413.528.7744 or 800.562.7569

thecgroup.com

VATS Trainer for Explants LEFT-SIDED

Accommodates porcine explants for skills development in thoracoscopic lobectomy and lung resection



Trainer Includes

- Rigid hemi-thorax in lateral decubitus position providing a shoulder landmark and access to the left anterior and posterior thoracic wall
- Incisable and replaceable skin/muscle element (#8013), allowing palpation of the ribs for port placement
- Hemi-ribcage with retractable ribs, derived from patient CT dataset
- Base with diaphragm, spine, and fluid drain
- Effluence tray
- Shipping case



*Bringing Practice to the
Practice of Medicine®*

thechamberlaingroup

INSTITUTIONAL DISCOUNTS AVAILABLE

TO ORDER CALL: +1.413.528.7744 or 800.562.7569 info@thecgroup.com

Quantity discounts may be available. Many of our products are made to order; please discuss any timing requirements with us in advance.
B-NP-171213 © 2017 The Chamberlain Group, LLC